



STARD10 Antibody

Product Code	CSB-PA022799GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y365
Immunogen	Human STARD10
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	StAR-related lipid transfer (START) domain containing 10;STARD10;CGI-52;MGC14401;NY-CO-28;PCTP2;SDCCAG28 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	STARD10
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.