



STARD8 Antibody

Product Code	CSB-PA022807GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92502
Immunogen	Human STARD8
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	StAR-related lipid transfer (START) domain containing 8;STARD8;DKFZp686H1668;DLC3;KIAA0189;STARTGAP3 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	STARD8
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.