



STK24 Antibody

Product Code	CSB-PA022841GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y6E0
Immunogen	Human STK24
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	serine/threonine kinase 24 (STE20 homolog, yeast);STK24;MST-3;MST3;MST3B;STE20;STK3 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	STK24
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.