



STK32C Antibody

Product Code	CSB-PA022847GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86UX6
Immunogen	Human STK32C
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	serine/threonine kinase 32C;STK32C;MGC23665;PKE;RP11-140A10.1;YANK3 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	STK32C
Target Details	This protein is a member of the serine/threonine protein kinase family. The specific function of this kinase is not known.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.