



# STK36 Antibody

|                            |  |
|----------------------------|--|
| <b>Product Code</b>        | CSB-PA022851GA01HU   |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                            |
| <b>Uniprot No.</b>         | Q9NRP7   |
| <b>Immunogen</b>           | Human STK36  |
| <b>Raised In</b>           | Rabbit   |
| <b>Species Reactivity</b>  | Human,Mouse,Rat  |
| <b>Tested Applications</b> | ELISA,WB,IHC   |
| <b>Storage Buffer</b>      | PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.     |
| <b>Purification Method</b> | Antigen Affinity Purified  |
| <b>Isotype</b>             | IgG  |
| <b>Alias</b>               | serine/threonine kinase 36, fused homolog (Drosophila);STK36;DKFZp434N0223;FU;KIAA1278 ; |
| <b>Product Type</b>        | Purified Rabbit Anti human PolyClonal Antibody   |
| <b>Immunogen Species</b>   | Homo sapiens (Human)   |
| <b>Target Names</b>        | STK36  |
| <b>Usage</b>               | For Research Use Only. Not for use in diagnostic or therapeutic procedures.              |