



STMN2 Antibody

Product Code	CSB-PA022858GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q93045
Immunogen	Human STMN2
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Form	Liquid
Storage Buffer	PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	STMN2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.