



STOML1 Antibody

Product Code	CSB-PA022863GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UBI4
Immunogen	Human STOML1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	stomatin (EPB72)-like 1;STOML1;FLJ36370;SLP-1;STORP;hUNC-24 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	STOML1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.