



STT3A Antibody

Product Code	CSB-PA022885GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P46977
Immunogen	Human STT3A
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	STT3, subunit of the oligosaccharyltransferase complex, homolog A (S. cerevisiae);STT3A;FLJ27038;ITM1;MGC9042;STT3-A;TMC ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	STT3A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.