



SUCLG2 Antibody

Product Code	CSB-PA022920GA01HU
Storage	PBS with 0.02% sodium azide and 50% glycerol, pH 7.3.
Uniprot No.	Q96I99
Immunogen	Human SUCLG2
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Form	Liquid
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SUCLG2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.