



SULT4A1 Antibody

Product Code	CSB-PA022944GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BR01
Immunogen	Human SULT4A1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	sulfotransferase family 4A, member 1;SULT4A1;BR-STL-1;BRSTL1;DJ388M5.3;MGC40032;NST;SULTX3;hBR-STL-1 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SULT4A1
Target Details	This gene encodes a member of the sulfotransferase family. The encoded protein is a brain-specific sulfotransferase believed to be involved in the metabolism of neurotransmitters. Polymorphisms in this gene may be associated with susceptibility to schizophrenia.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.