



SUV39H2 Antibody

Product Code	CSB-PA022973GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H511
Immunogen	Human SUV39H2
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	suppressor of variegation 3-9 homolog 2 (Drosophila);SUV39H2;FLJ23414;KMT1B ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SUV39H2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.