



SV2B Antibody

Product Code	CSB-PA022979GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7L1I2
Immunogen	Human SV2B
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	synaptic vesicle glycoprotein 2B;SV2B;HsT19680;KIAA0735 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SV2B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.