



SYCE1 Antibody

Product Code	CSB-PA022987GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N0S2
Immunogen	Human SYCE1
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	synaptonemal complex central element protein 1;SYCE1;C10orf94;CT76;RP11-108K14.6;bA108K14.6 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SYCE1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.