



SYNCRIP Antibody

Product Code	CSB-PA023007GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60506
Immunogen	Human SYNCRIP
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	synaptotagmin binding, cytoplasmic RNA interacting protein;SYNCRIP;GRY-RBP;HNRPQ1;NSAP1;RP1-3J17.2;dJ3J17.2;hnRNP-Q;pp68 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	SYNCRIP
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.