



SYTL4 Antibody

Product Code	CSB-PA023049GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96C24
Immunogen	Human SYTL4
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	synaptotagmin-like 4;SYTL4;DKFZp451P0116;FLJ40960;SLP4 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SYTL4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.