



TAX1BP3 Antibody

Product Code	CSB-PA023180GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14907
Immunogen	Human TAX1BP3
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	Tax1 (human T-cell leukemia virus type I) binding protein 3;TAX1BP3;TIP-1 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TAX1BP3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.