



TCF3 Antibody

Product Code	CSB-PA023301GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P15923
Immunogen	Human TCF3
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	transcription factor 3 (E2A immunoglobulin enhancer binding factors E12/E47);TCF3;E2A;ITF1;MGC129647;MGC129648;bHLHb21 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TCF3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.