



TCP11 Antibody

Product Code	CSB-PA023324GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WWU5
Immunogen	Human TCP11
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	t-complex 11 homolog (mouse);TCP11;D6S230E;KIAA0229;MGC111103 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TCP11
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.