



TCTN2 Antibody

Product Code	CSB-PA023339GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96GX1
Immunogen	Human TCTN2
Species Reactivity	Human, Mouse, Rat
Tested Applications	ELISA, IF, IP
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	697aa, 77kd;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TCTN2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.