



TEX13A Antibody

Product Code	CSB-PA023403GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BXU3
Immunogen	Human TEX13A
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	testis expressed 13A; TEX13A; MGC131984 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TEX13A
Target Details	This gene is similar to a mouse gene that is expressed in the testis.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.