



THNSL2 Antibody

Product Code	CSB-PA023500GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86YJ6
Immunogen	Human THNSL2
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	threonine synthase-like 2 (<i>S. cerevisiae</i>);THNSL2;FLJ10916;THS2;TSH2 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	THNSL2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.