



TICAM2 Antibody

Product Code	CSB-PA023533GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86XR7
Immunogen	Human TICAM2
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	toll-like receptor adaptor molecule 2; TICAM2; MGC129876; MGC129877; TICAM-2; TIRAP3; TIRP; TRAM ;
Product Type	Purified Rabbit Anti Human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TICAM2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.