



TIMM13 Antibody

Product Code	CSB-PA023549GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y5L4
Immunogen	Human TIMM13
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	translocase of inner mitochondrial membrane 13 homolog (yeast);TIMM13;TIM13;TIM13B;TIMM13A;TIMM13B;ppv1 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TIMM13
Target Details	This gene encodes a translocase with similarity to yeast mitochondrial proteins that are involved in the import of metabolite transporters from the cytoplasm and into the mitochondrial inner membrane. The encoded protein and the TIMM8a protein form a 70 kDa complex in the intermembrane space. This gene is in a head-to-tail orientation with the gene for lamin B2.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.