



# TIMM13 Antibody

|                            |                                                                                                                                                                                                                                                                                                                                                                              |
|----------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>        | CSB-PA023549GA01HU                                                                                                                                                                                                                                                                                                                                                           |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                                                                                                                                                                                                                                                                                                                |
| <b>Uniprot No.</b>         | Q9Y5L4                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Immunogen</b>           | Human TIMM13                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Raised In</b>           | Rabbit                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Species Reactivity</b>  | Human,Mouse,Rat                                                                                                                                                                                                                                                                                                                                                              |
| <b>Tested Applications</b> | ELISA,WB                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Storage Buffer</b>      | PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.                                                                                                                                                                                                                                                                                         |
| <b>Purification Method</b> | Antigen Affinity Purified                                                                                                                                                                                                                                                                                                                                                    |
| <b>Isotype</b>             | IgG                                                                                                                                                                                                                                                                                                                                                                          |
| <b>Alias</b>               | translocase of inner mitochondrial membrane 13 homolog (yeast);TIMM13;TIM13;TIM13B;TIMM13A;TIMM13B;ppv1 ;                                                                                                                                                                                                                                                                    |
| <b>Product Type</b>        | Purified Rabbit Anti human PolyClonal Antibody                                                                                                                                                                                                                                                                                                                               |
| <b>Immunogen Species</b>   | Homo sapiens (Human)                                                                                                                                                                                                                                                                                                                                                         |
| <b>Target Names</b>        | TIMM13                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Target Details</b>      | This gene encodes a translocase with similarity to yeast mitochondrial proteins that are involved in the import of metabolite transporters from the cytoplasm and into the mitochondrial inner membrane. The encoded protein and the TIMM8a protein form a 70 kDa complex in the intermembrane space. This gene is in a head-to-tail orientation with the gene for lamin B2. |