



TLE3 Antibody

Product Code	CSB-PA023587GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q04726
Immunogen	Human TLE3
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	transducin-like enhancer of split 3 (E(sp1) homolog, Drosophila);TLE3;ESG;ESG3;FLJ39460;GRG3;HsT18976;KIAA1547 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TLE3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.