



TMEM173 Antibody

Product Code	CSB-PA023754GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86WV6
Immunogen	Human TMEM173
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC, IF, FC
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	transmembrane protein 173; TMEM173; FLJ38577; MITA; MPYS; STING ;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TMEM173
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.