



TMEM50B Antibody

Product Code	CSB-PA023850GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P56557
Immunogen	Human TMEM50B
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	transmembrane protein 50B;TMEM50B;C21orf4;DKFZp686C2482;FLJ26146;HCVP7TP3 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TMEM50B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.