



TMPO Antibody

Product Code	CSB-PA023913GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P42166
Immunogen	Human TMPO
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC,IF
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	thymopoietin;TMPO;CMD1T;LAP2;LEMD4;MGC61508;PRO0868;TP ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TMPO
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.