



TMX2 Antibody

Product Code	CSB-PA023951GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y320
Immunogen	Human TMX2
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	thioredoxin-related transmembrane protein 2;TMX2;CGI-31;DKFZp781O2021;MGC111151;PDIA12;PIG26;TXNDC14 ;
Product Type	Purified Rabbit Anti Human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TMX2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.