



TOM1L1 Antibody

Product Code	CSB-PA024044GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75674
Immunogen	Human TOM1L1
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	target of myb1 (chicken)-like 1; TOM1L1; OK/KNS-CL.3; SRCASM ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TOM1L1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.