



TPPP Antibody

Product Code	CSB-PA024115GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O94811
Immunogen	Human TPPP
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	tubulin polymerization promoting protein; TPPP; TPPP/p25; TPPP1; p24; p25; p25alpha ;
Product Type	Purified Rabbit Anti Human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TPPP
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.