



TRAPPC9 Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA024234GA01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q96Q05 |
| Immunogen | Human TRAPPC9 |
| Raised In | Rabbit |
| Species Reactivity | Human,Mouse,Rat |
| Tested Applications | ELISA,WB,IHC |
| Storage Buffer | PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles. |
| Purification Method | Antigen Affinity purified |
| Isotype | IgG |
| Alias | trafficking protein particle complex 9;TRAPPC9;IBP;IKBKBBP;KIAA1882;MGC4737;MGC4769;NIBP;T1;TRS120 ; |
| Product Type | Purified Rabbit Anti Human,Mouse PolyClonal Antibody |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | TRAPPC9 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |