



Recombinant *Saccharomyces cerevisiae* Cytoplasmic tRNA 2-thiolation protein 1 (NCS6)

Product Code	CSB-EP002426STA-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A6ZTX8
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain YJM789) (Baker's yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MSFTAPSDPV NKPTKVKVSQ LCELCHSRKA LIRRPKNLSK LCKQCFLVF ETEIHNTIVA NNLFQRGEKV AVGASGGKDS TVLAHMLKLL NDRYDYGIEI VLLSIDEGLI GYRDDSLATV KRNNQQYGLP LEIFSFKDLY DWTMDEIVSV AGIRNSCTYC GVFRRLQSLDR GAAKLGISHV VTGHNADDMA ETVLMNILRG DVARLEKSTA IITQSSGSPV KRSKPFKYSY QKEIVLYAHY MKLDYFSTEC TYAPEAFRGT AREYMKNLEA VRPSCIIDII QSGENLALKA KKSNAKRKRVV KFVDGNRCAR CGYLSSNNIC KACMLLEGLE KSRAQVAIEN DTSADGAALK LRALEKLSF
Source	<i>E.coli</i>
Target Names	NCS6
Protein Names	Recommended name: Cytoplasmic tRNA 2-thiolation protein 1 EC= 2.7.7.- Alternative name(s): Cytoplasmic tRNA adenylyltransferase 1 Needs CLA4 to survive protein 6 Thiolation of uridine in cytoplasmic tRNA protein 1
Expression Region	1-359
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.