



TRIM46 Antibody

Product Code	CSB-PA024484GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7Z4K8
Immunogen	Human TRIM46
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IF
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	tripartite motif-containing 46;TRIM46;FLJ23229;GENEY;TRIFIC ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TRIM46
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.