



TRPC1 Antibody

Product Code	CSB-PA025075GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P48995
Immunogen	Human TRPC1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	transient receptor potential cation channel, subfamily C, member 1;TRPC1;HTRP-1;MGC133334;MGC133335;TRP1 ;
Product Type	Purified Rabbit Anti Human,Mouse PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TRPC1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.