



TRPC4AP Antibody

Product Code	CSB-PA025079GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TEL6
Immunogen	Human TRPC4AP
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	transient receptor potential cation channel, subfamily C, member 4 associated protein;TRPC4AP;C20orf188;TRRP4AP;TRUSS ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TRPC4AP
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.