



TSC22D4 Antibody

Product Code	CSB-PA025118GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y3Q8
Immunogen	Human TSC22D4
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	0.1M NaHCO ₃ ,0.1M Glycine, 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	TSC22 domain family, member 4;TSC22D4;THG-1;THG1 ;
Product Type	Purified Rabbit Anti Human,Mouse PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TSC22D4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.