



TXK Antibody

Product Code	CSB-PA025362GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P42681
Immunogen	Human TXK
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,IHC,IF
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	TXK tyrosine kinase;TXK;BTKL;MGC22473;PSCTK5;PTK4;RLK;TKL ;
Product Type	Purified Rabbit Anti Human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TXK
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.