



UBE2O Antibody

Product Code	CSB-PA025473GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9C0C9
Immunogen	Human UBE2O
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	ubiquitin-conjugating enzyme E2O;UBE2O;E2-230K;FLJ12878;KIAA1734 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	UBE2O
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.