



UBL7 Antibody

Product Code	CSB-PA025503GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96S82
Immunogen	Human UBL7
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	ubiquitin-like 7 (bone marrow stromal cell-derived);UBL7;BMSC-UbP;MGC14421;TCBA1 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	UBL7
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.