



UBR7 Antibody

Product Code	CSB-PA025521GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N806
Immunogen	Human UBR7
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	ubiquitin protein ligase E3 component n-recognin 7 (putative);UBR7;C14orf130;FLJ10483;MGC9518 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	UBR7
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.