



UCHL5 Antibody

Product Code	CSB-PA025544GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y5K5
Immunogen	Human UCHL5
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	ubiquitin carboxyl-terminal hydrolase L5;UCHL5;CGI-70;INO80R;UCH37 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	UCHL5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.