



USP34 Antibody

Product Code	CSB-PA025725GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q70CQ2
Immunogen	Human USP34
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	ubiquitin specific peptidase 34; USP34; FLJ43910; KIAA0570; KIAA0729; MGC104459 ;
Product Type	Purified Rabbit Anti Human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	USP34
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.