



VPS37A Antibody

Product Code	CSB-PA025906GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NEZ2
Immunogen	Human VPS37A
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	vacuolar protein sorting 37 homolog A (S. cerevisiae);VPS37A;FLJ32642;FLJ42616;HCRP1;PQBP2 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	VPS37A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.