



VPS53 Antibody

Product Code	CSB-PA025917GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5VIR6
Immunogen	Human VPS53
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	vacuolar protein sorting 53 homolog (S. cerevisiae);VPS53;FLJ10979;HCCS1;MGC39512;hVps53L;pp13624 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	VPS53
Target Details	This gene encodes a protein with sequence similarity to the yeast Vps53p protein. Vps53p is involved in retrograde vesicle trafficking in late Golgi.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.