



VTA1 Antibody

Product Code	CSB-PA025940GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NP79
Immunogen	Human VTA1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IF
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	Vps20-associated 1 homolog (S. cerevisiae);VTA1;C6orf55;DRG-1;DRG1;FLJ27228;HSPC228;LIP5;My012;SBP 1 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	VTA1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.