



WBSCR16 Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA025981GA01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q96I51 |
| Immunogen | Human WBSCR16 |
| Raised In | Rabbit |
| Species Reactivity | Human,Mouse,Rat |
| Tested Applications | ELISA,WB,IHC |
| Storage Buffer | PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles. |
| Purification Method | Antigen Affinity Purified |
| Isotype | IgG |
| Alias | Williams-Beuren syndrome chromosome region 16;WBSCR16;DKFZp434D0421;MGC189739;MGC44931 ; |
| Product Type | Purified Rabbit Anti human PolyClonal Antibody |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | RCC1L |
| Target Details | This gene encodes an RCC1-like G-exchanging factor. It is deleted in Williams syndrome, a multisystem developmental disorder caused by the deletion of contiguous genes at 7q11.23. |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |