



WIPI2 Antibody

Product Code	CSB-PA026118GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y4P8
Immunogen	Human WIPI2
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	WD repeat domain, phosphoinositide interacting 2;WIPI2;Atg21;CGI-50;DKFZp434J154;DKFZp686P02188;FLJ12979;FLJ14217;FLJ42984;WIPI-2 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	WIPI2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.