



# Recombinant Human Zinc finger protein 143 (ZNF143)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-MP026542HU  |
| <b>Storage</b>           | Store at -20°C, for extended storage, conserve at -20°C or -80°C.   |
| <b>Uniprot No.</b>       | P52747  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | Homo sapiens (Human)  |
| <b>Purity</b>            | ≥85% (SDS-PAGE)   |
| <b>Sequence</b>          | MLLAQINRDS QGMTEFPGGG MEAQHVTLCL TEAVTVADGD NLENMEGVSL<br>QAVTLADGST AYIQHNSKDA KLIDGQVIQL EDGSAAYVQH VPIPKSTGDS<br>LRLEDGQAVQ LEDGTTAFIH HTSKDSYDQS ALQAVQLEDG TTAYIHHAVQ<br>VPQSDTILAI QADGTVAGLH TGDATIDPDT ISALEQYAAK VSIDGSESVA<br>GTGMIGENEQ EKKMQIVLQG HATRVTAKSQ QSGEKAFCRCE YDGCCKLYTT<br>AHLKVVHERS HTGDRPYQCE HAGCGKAFAT GYGLKSHVRT HTGEKPYRCS<br>EDNCTKSFKT SGDLQKHIRT HTGERPFKCP FEGCGRSFTT SNIRKVHVRT<br>HTGERPYCT EPGCGRAFAS ATNYKNHVRI HTGEKPYVCT VPGCDKRFTE<br>YSSLYKHHVV HTHSKPYN CN HCGKTYKQIS TLAMHKRTAH NDTEPIEEEEQ<br>EAFFEP PPQ GEDVLKGSQI TYVTGVEGDD VVSTQVATVT QSGLSQQVTL<br>ISQDGTQHVN ISQADMQAIG NTITMVTQDG TPITVPAHDA VISSAGTHSV<br>AMVTAEGTEG EQVAIVAQDL AAFHTASSEM GHQQHSHHLV TTETRPLTLV<br>ATSNGTQIAV QLGEQPSLEE AIRIASRIQQ GETPGLDD |
| <b>Source</b>            | Mammalian cell  |
| <b>Target Names</b>      | ZNF143  |
| <b>Protein Names</b>     | Recommended name: Zinc finger protein 143 Alternative name(s): SPH-binding factor Selenocysteine tRNA gene transcription-activating factor Short name= hStaf  |
| <b>Expression Region</b> | 1-638   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | Full length protein   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.   |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life  |



of lyophilized form is 12 months at -20°C/-80°C.