



YBX1 Antibody

Product Code	CSB-PA026247GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P67809
Immunogen	Human YBX1
Raised In	Rabbit
Species Reactivity	Human,Mouse
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	Y box binding protein 1;YBX1;BP-8;CSDA2;CSDB;DBPB;MDR-NF1;MGC104858;MGC110976;MGC117250;NSEP-1;NSEP1;YB-1;YB1 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	YBX1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.